
AFDM411 Electric Controlled Continuous Zoom and Autofocus Digital Microscope



Contents

1	Introduction to AFDM411	1
1.1	The Module Specifications Of AFDM411	1
1.1.1	<i>AFDM Camera Module Datasheet.....</i>	<i>2</i>
1.1.2	<i>AFDM Lens Module Datasheet</i>	<i>2</i>
1.1.3	<i>AFDM Light Module.....</i>	<i>2</i>
1.2	AFDM411 Characteristic And Specification	2
1.2.1	<i>The Basic Characteristic of AFDM411.....</i>	<i>2</i>
1.2.2	<i>Specification of AFDM411.....</i>	<i>3</i>
1.2.3	<i>Dimension of AFDM411</i>	<i>5</i>
1.2.4	<i>Packing Information of AFDM411</i>	<i>5</i>
2	Installation and Operation of AFDM Series Product	6
3	Images Captured with AFDM411	7
4	Software and App	10
5	AFDM411 Series Camera Configurations	11
5.1	Camera Working Standalone with Built-in XFCAMView Software.....	11
5.2	Connecting Camera to Computers with USB2.0 Port	12
5.3	Camera Working in WLAN Mode (AP Mode).....	13
5.4	Connecting Camera to The PC with LAN Port.....	15
5.5	Connecting Multi-Cameras to The Router Through The LAN Port/WLAN STA Mode for The Network Application	16
6	Introduction of XFCAMView UI and Functions.....	20
6.1	Control UI.....	20
6.2	The Camera Control Panel on the Left Side of the Video Window	21
6.3	The Measurement Toolbar On The Upper Side Of The Video Window	21
6.3.1	<i>Introduction to Measurement Toolbar</i>	<i>21</i>
6.3.2	<i>Calibration Method</i>	<i>22</i>
6.4	Synthesis Camera Control Toolbar At The Bottom Of The Video Window	23
6.4.1	<i>Setting>Network.....</i>	<i>24</i>
6.4.2	<i>Setting>Measurement.....</i>	<i>26</i>
6.4.3	<i>Setting>Magnification.....</i>	<i>26</i>
6.4.4	<i>Settings>Image Format</i>	<i>27</i>
6.4.5	<i>Setting>Video.....</i>	<i>27</i>
6.4.6	<i>Setting>Storage</i>	<i>28</i>
6.4.7	<i>Setting>Files.....</i>	<i>28</i>
6.4.8	<i>Setting>Time.....</i>	<i>29</i>
6.4.9	<i>Setting>Language.....</i>	<i>29</i>
6.4.10	<i>Setting>Miscellaneous</i>	<i>30</i>
6.4.11	<i>Auto Focus Control Panel on the Right Side of Video Window</i>	<i>30</i>
6.4.12	<i>Focus Region On The Video Window</i>	<i>31</i>

1 Introduction to AFDM411

AFDM is a series of electric controlled continuous zoom and autofocus all-in-one digital microscope with a large field of view by Touptek Photonics. It is integrated with HDMI camera, Electric Controlled Continuous Zoom Auto-focus Objective and LED Integrated Illumination Light. AFDM is the abbreviation of Auto-focus Digital Microscope. Different products in the AFDM series can be formed with different part to satisfy the application requirement.

AFDM can be assembled with various brackets or arms and offer a continuous zooming ratio with different lens. AFDM also supports autofocus mode and manual focus mode.

AFDM comes with a high-performance SONY CMOS sensor. It also has an embedded ARM core, allowing the camera to be connected directly to the HDMI monitor. The camera has XFCAMView software built within it, including Camera Control Panel, Auto Focus Control Panel, Measurement Toolbar, and Synthesis Camera Control Toolbar. Users can directly control the camera and perform various operations through a USB mouse. The images and videos captured by AFDM can be saved on an SD card for on-site analysis and follow-up research.

AFDM can be widely used in industrial inspection, medical observation, teaching and scientific research, automation system, and other fields.

AFDM101 supports 1080P60FPS HDMI output.

AFDM411 supports 4K30FPS HDMI output and USB/ETH/WIFI outputs.



Figure 1-1 AFDM's Front and Back View



Figure 1-2 AFDM's Side and Front(with LED light) View

1.1 The Module Specifications Of AFDM411

1.1.1 AFDM Camera Module Datasheet

Order Code	Sensor & Size(mm)	Pixel(μm)	G Sensitivity/Dark Signal	FPS/Resolution	Binning	Exposure(ms)
HI080PA	Sony IMX462LQR-C 1/2.8"(5.57x3.13)	2.9x2.9	921mv/0.15mv with 1/30s	60/1920*1080(HDMI)	1x1	0.06~918
H4KPA	Sony IMX415LQR-C 1/2.8"(5.57x3.13)	1.45x1.45	300mv/0.13 with 1/30s	30@3840*2160(HDMI) 30@3840*2160(NETWORK) 30@3840*2160(USB)	1x1	0.04~1000

C: Color; M: Monochrome;

1.1.2 AFDM Lens Module Datasheet

Order Code	Working Distance(mm)	Zoom Range	MTF(lp/mm)	Distortion	FOV@1X(mm)	FOV@20X(mm)
EMZO-20XA	150~195	0.028X~0.56X	160	0.5%	200x112.5	10x5.6

1X and 20x are defined as the normalized magnification, which is only used to represent the relative relationship between the lowest and highest magnification. Here, the normalized equations are $1x = 0.028/0.028$; $20X=0.56/0.028$;

1.1.3 AFDM Light Module

Order Code	LED	Power	Inner Dia.(mm)	Out Dia.(mm)		
DRL-5076A-NPC	8 CREE xpes	3V/3A	50	76		

DRL: LED direct ring light with adjustable brightness; NPC: No power cable

1.2 AFDM411 Characteristic And Specification

The [AFDM411](#) comes with [H4KPA HDMI](#) camera, [EMZO-20XA](#) lens and [DRL-5076A-NPC](#) light source(Optional);

1.2.1 The Basic Characteristic of AFDM411

- 5 groups 16 elements EMZO with 0.028~0.56X, 20 zoom ratio, supports auto and manual focus
- 192mm standard working distance with 150~195mm depth of field
- At standard working distance, the large field of view 200mm*112.5mm at low magnification, helping users to quickly locate the target object, the small field of view 10mm*5.6mm at higher magnification, helping users to observe microscopically
- Sony 1/2.8" 4K Starvis CMOS with high signal-to-noise ratio
- 4K HDMI/USB/ETH/WiFi multiple video outputs
- 4K/1080P auto switching according to monitor resolution
- SD card/USB flash drive for captured image and video storage, support local preview and playback
- Built-in mouse control software [XFCAMView](#), all functions can be realized with USB mouse
- Embedded mouse Camera Control Panel, Measurement Toolbar, Synthesis Control Toolbar, AF Control Panel
- Multi-language support
- Head suction LED ring light, the brightness can be directly controlled by [XFCAMView](#)
- With the adapter bracket of 76mm diameter, a electric controlled continuous zoom AFDM can be built



Figure 1-3 TPS-30A(bracket)+AFDM411+1080P Monitor

1.2.2 Specification of AFDM411



Interface & Button Functions	
USB Mouse	USB mouse for XFCAMView control
USB2.0	Connect USB flash drive to save pictures and videos Connect 5G WLAN module to transfer video wirelessly in real time with ToupView/ToupLite
HDMI	Comply with HDMI1.4 standard. 4K/1080P format video output and supporting automatic switch between 4K and 1080P format according to the connected monitors
USB Video	Connect PC or other host device to realize video image transmission with ToupView/ToupLite
LAN	LAN port to connect router and switch to transfer video with ToupView/ToupLite
ON/OFF	Power on/off switch
LED	Power LED indicator
SD	Comply with SDIO3.0 standard and SD card could be inserted for video and images saving
DC12V3A	DC12V3A power input
XFCAMView Software Functions	
UI Operation	With USB mouse to operate on the embedded XFCAMView

Image Capture	8M (3840*2160) JPEG/TIFF image in SD card or USB flash drive
Video Record	Video format: 8M(3840*2160) H264/H265 encoded MP4 file Video saving frame rate:30fps
Camera Control Panel	Including Exposure, Gain, White Balance, Sharpness, Denoise, Denoise, Saturation, Gamma, Contrast, Brightness, Power Frequency control
Measurement Toolbar	Including Calibration, Measurement , and measurement parameter Export functions
Synthesis Control Toolbar	Including software Zoom, Flip, Freeze, Crosshair, LED Control, Auto-focus, Comparison, Browser, Setting, Version Check function
Auto Focus Control Panel	Including Zoom, Auto Focus, One Push, Manual Focus, Reset , and other functions
Software ToupView/ToupLite Environment under LAN/WLAN/USB Video Output	
White Balance	Auto White Balance
Color Technique	Ultra-Fine Color Engine
Capture/Control SDK	Windows/Linux/macOS/Android Multiple Platform SDK(Native C/C++, C#/VB.NET, Python, Java, DirectShow, Twain, etc)
Recording System	Still Picture or Movie
Operating System	Microsoft Windows XP / Vista / 7 / 8 / 8.1 /10(32 & 64 bit)/ ToupView OSx(Mac OS X)/ ToupLite Linux/ ToupLite
PC Requirements	CPU: Equal to Intel Core2 2.8GHz or Higher
	Memory: 4GB or More
	Ethernet Port: RJ45 Ethernet Port
	Display:19" or Larger
	CD-ROM
Operating Environment	
Operating Temperature(in Centidegree)	-10~ 50
Storage Temperature(in Centidegree)	-20~ 60
Operating Humidity	30~80%RH
Storage Humidity	10~60%RH
Dimension	
Length x Width x Height	80mm x 80mm x 80mm
Shipping Weight	0.75kg

1.2.3 Dimension of AFDM411



Figure 1-4 Dimension of AFDM411

1.2.4 Packing Information of AFDM411



Figure 1-5 Packing Information of AFDM411

Standard Packing List		
A	Gift box: L:220cm W:220cm H:110cm (1pcs, 2.0kg/box)	
B	AFDM411	
C	Power Adapter: Input: AC 100~240V 50Hz/60Hz, Output: DC 12V 3A	American Standard: Model: HKA03612030-7K : UL/CE/FCC(With American Standard AC Power Cable) European Standard: Model: HKA03612030-7K : UL/CE/FCC(With European Standard AC Power Cable) EMI Standard: FCC Part 15 Subpart B EMS Standard: EN61000-4-2,3,4,5,6


D	USB Mouse	
E	HDMI Cable	
F	USB2.0 A male to A male gold-plated connectors cable /2.0m	
G	CD (Driver & utilities software, Ø12cm)	
Optional Accessory		
H	Ethernet cable	
I	LED Ring Light(DRL-5076A-NPC)	
J	USB flash drive	
K	USB WLAN adapter	
L	SD card(16G)	
M	Calibration kit	106011/TS-M1(X=0.01mm/100Div.); 106012/TS-M2(X,Y=0.01mm/100Div.); 106013/TS-M7(X=0.01mm/100Div., 0.10mm/100Div.)

2 Installation and Operation of AFDM Series Product

Before use, please install the [AFDM](#) series product on an adaptive bracket.

- 1.Plug HDMI cable into the [HDMI](#) port to connect [AFDM](#) and HDMI monitor;
- 2.Plug a USB mouse into [USB Mouse](#) port, to get control of the [AFDM](#) by using built-in software [XFCAMView](#);
- 3.Plug DC12V3A power adapter into [DC12V3A](#) port, to supply power for the [AFDM](#), the [LED Indicator](#) will turn into red;
- 4.Insert SD card into [SD card Slot](#) for saving captured images and recorded videos;
- 5.Press [ON/OFF](#) button to start the [AFDM](#), [LED Indicator](#) will turn into blue;
- 6.Move mouse to the left side of the video window, the [Camera Control Panel](#) will appear. It includes [Manual/Automatic Exposure](#), [White Balance](#), [Sharpness](#), [Denoise](#), and other functions, please refer to section 3.2 for details;
- 7.Move mouse to the upper side of the video window, the [Measurement Toolbar](#) will appear. It includes [calibration](#), measurement of [lines](#), [angles](#), [rectangles](#), [circles](#), etc, and supports data export(*.CSV format), please refer to section 3.3 for detail;
- 8.Move mouse to the bottom side of the video window, the [Synthesis Camera Control Toolbar](#) will appear. Operations like [Zoom In](#), [Zoom Out](#), [Flip](#), [Freeze](#), [Crossline](#), [LED brightness control](#), [Autofocus](#), [SD card contents browsing](#), [Settings](#), and [Camera Version](#) can be executed. See section 3.4 for details;
- 9.Move mouse to the bottom side of the video window, the [Synthesis Camera Control Toolbar](#) will pop up automatically. Click [AF](#) button, and [Auto Focus Control Panel](#) will show up for autofocus operation, it supports 20X optical zoom, [Autofocus](#), [Manual Focus](#), [Reset](#), and [One Push](#) operation. See section 3.5 for details.


3 Images Captured with AFDM411



Ruler Captured with AFDM411 at 1X




Ruler Captured with AFDM411 at 10X




Ruler Captured with AFDM411 at 20X




Print Captured with AFDM411 at 1.0X




Print Captured with AFDM411 at 10X




Print Captured with AFDM411 at 20X



PCB Captured with AFDM411 at 4.0X



PCB Captured with AFDM411 at 10X



PCB Captured with AFDM411 at 20X

4 Software and App

The software or the **APP** can be downloaded from the following link:

Windows: <https://www.touptekphotonics.com/download/>

Linux & macOS: <https://www.touptekphotonics.com/download/>

iOS: <https://itunes.apple.com/us/app/toupview/id911644970>

Android: <https://play.google.com/store/apps/details?id=com.touptek.tpview>

For **ToupLite** and **ToupView App**, the **Auto-focus** and **LED Brightness Control** are not available


5 AFDM411 Series Camera Configurations

You can use the AFDM411 series camera in 5 different ways. Each connection requires different hardware configuration.


5.1 Camera Working Standalone with Built-in XFCAMView Software

For this application, apart from the microscope, you only need an HDMI monitor, the supplied USB mouse, and the camera embedded with XFCAMView software. The steps to start the camera are listed as below:


- Connect the camera to a HDMI monitor using the HDMI cable;




- Insert the supplied USB mouse to the camera's USB port;



- Insert the supplied SD card/USB flash drive (USB2.0 slot) into the AFDM411 series camera SD card slot/USB2.0 slot;



- Connect the camera to the power adapter and turn it on;



- Turn on the monitor and view the video in the [XFCAMView](#) software. Move the mouse to the left, top or bottom of the [XFCAMView](#) UI, different control panel or toolbar will pop up and users could operate with the mouse at ease.





Figure 5-1 XFCAMView and AFDM411 Series Camera in HDMI Mode

5.2 Connecting Camera to Computers with USB2.0 Port


For Windows user (Windows XP (32bit), Windows 7/8/10/11 (32/64 bit)), please use [ToupView](#).

For [macOS](#) and [Linux](#) user (macOS 10.10 or above or [Linux](#) distributions with kernel 2.6.27 or higher), please use [ToupLite](#). The steps to start the camera are listed below:


- Start the camera according to Sec. 5.1. After the camera is running, connect camera to computer with USB cable. Please use “USB Video” slot, not “USB Mouse” slot as shown below.
- Install [ToupView/ToupLite](#) on your PC or install [ToupView App](#) on the mobile device; Run the software [ToupView/ToupLite](#), clicking the camera name in the camera list n to start the live video as shown in Figure 5-2.



Plug the [USB WLAN](#) adapter into the camera's USB2 .0 port;



Install [ToupView/ToupLite](#) on your PC or install [ToupView App](#) on the mobile device, Connect the PC or mobile device to the camera's [WLAN AP](#) point; The network name (SSID) and the [WLAN](#) password (The default one is 12345678) can be found on the camera's [Setting>Network>WLAN](#) page in [AP](#) mode.



Start [ToupView/ToupLite](#) software or [ToupView App](#) and check the configuration. Normally, the active AFDM411 series cameras will be automatically recognized. The live image of each camera is shown in Figure 5-3. For the display, the [Camera List](#) tool window is used in [ToupView/ToupLite](#) software, and the [Camera Thumbnail](#) is used in [ToupView App](#).





Figure 5-3 ToupView and AFDM411 Series Camera in WLAN AP Mode

5.4 Connecting Camera to The PC with LAN Port

This application uses the camera as the network camera. User must configure the IP of the camera and PC manually and ensure their IP addresses are in the same net. The subnet mask and gateway of the camera and PC must be the same.

Start the camera according to Sec. 5.1 after the camera is running, clicking  button on the [Synthesis Camera Control Toolbar](#) at the bottom of the video window (See Figure 5-1), a small window called [Settings](#) will pop up as shown below on the left side, clicking [LAN](#) property page, uncheck the DHCP item. Input [IP Address](#), [Subnet Mask](#) and [Default Gateway](#) for the camera. Designate [Internet Protocol Version 4 \(TCP/IPv4\) Settings](#) page's IP address on the PC with similar configuration as shown below on the right side but with different [IP address](#).




Figure 5-4 Configure the AFDM411 Series Camera IP






Figure 5-5 Configure the PC's IP

After the above configurations are finished, user can connect the AFDM411 series camera to the computer through the Ethernet cable as shown below:

Connect the LAN port with the Ethernet cable to the PC's network port;




Insert the supplied SD card/USB flash drive (USB2.0 slot) into the AFDM411 series camera's SD card slot/USB2.0 slot;



Install [ToupView/ToupLite](#) on your PC or install [ToupView App](#) on the mobile device; Run the software [ToupView/ToupLite](#), clicking the camera name in the camera list starts the live video as shown in Figure 5-3.

5.5 Connecting Multi-Cameras to The Router Through The LAN Port/WLAN STA Mode for The Network Application


In LAN/WLAN STA mode, the camera connects to the router by LAN port/WLAN STA mode. If a router with LAN/WLAN capability is used, users could connect the router with Ethernet cable/WLAN to control the camera.




The connection and configuration are just the same as in Sec.5.1 or Sec. 5.4. But here, users need to check [DHCP](#). If [Multicast](#) is disabled or is not supported, users should only select [Unicast](#). If [Multicast](#) is supported by the network, users could select [Multicast](#) to achieve a better performance, especially in the case that multi-users connecting to the same camera.

In addition, please guarantee that the broadcasting function is enabled in the network.


Active AFDM411 series camera is recognized by [ToupView/ToupLite](#) software or [ToupView App](#) and they are displayed as a camera list or thumbnail in the software or app as shown in Figure 5-2.

Or start the camera according to Sec. 5.1. After the camera is running, move the mouse to the bottom of the video window and clicking the  button on the [Synthesis Camera Control Toolbar](#) at the bottom of the video window, a small window called [Settings](#) will pop up as shown below. Clicking [Network>WLAN](#) property page and choosing the [STA](#) in the [Wi-Fi Mode](#) edit box(The factory default configuration is [AP mode](#)). Input the to be connected router's [SSID](#) and [Password](#) as shown below:




Install [ToupView /ToupLite](#) software on your PC. Alternatively, install the free [ToupView App](#) on the mobile device;


Plug the Ethernet cable into the camera's [LAN](#) port and the other end to the PC (for those connected to router with [WLAN STA](#) mode);



Or plug the [USB WLAN](#) adapter into the camera's USB2.0 port(for those connected to router with [WLAN STA](#) mode);



Finally, as shown below, 2 AFDM411 series cameras are connected to the router with LAN cable and 2 AFDM411 series cameras are connected to the same router with WLAN STA mode(The number of the cameras, the connection mode(LAN or WLAN STA)) connected to the router are determined by the router performance)



Make sure that your PC or your mobile device is connected to the LAN or WLAN of the router; Start [ToupView/ToupLite](#) software or [ToupView App](#) and check the configuration. Normally, active AFDM411 series cameras are automatically recognized. The live image of each camera is displayed. For the display, [Camera List](#) control panel window is used in [ToupView/ToupLite](#) software, and [Camera Thumbnail](#) is used in [ToupView App](#); Select the AFDM411 series camera you are interested in. To do so, double click the camera's name in [Camera List](#) tool window if you use [ToupView /ToupLite](#) software; If you use [ToupView App](#), tap the camera's thumbnail in [Camera List](#) page(See Figure 5-6)

About the routers/switches

It is suggested that routers/switches supporting 802.11ac 5G segment should be selected to achieve better wireless connection experience.

6 Introduction of XFCAMView UI and Functions

6.1 Control UI










Figure 6-1 XFCAMView and Its Control UI

AFDM's XFCAMView software operation UI is shown in Figure 6-1. It includes **Camera Control Panel** on the left side of the video window, **Measurement Toolbar** on the top of the video window, **Synthesis Camera Control Toolbar**, and **Autofocus Control Panel** on the right side of the video window.

Software Toolbar / Control Bar / Control Panel	
1	Move the mouse to the left side of the video window, the Camera Control Panel will pop up automatically;
2	Move the mouse to the bottom of the video window, the Synthesis Camera Control Toolbar will pop up automatically;
3	Move the mouse to the bottom of the video window, the Synthesis Camera Control Toolbar will pop up automatically. Click the  button and the Auto Focus Control Panel will appear for autofocus operation;
4	Move the mouse to the upper side of the video window, the Measurement Toolbar will pop up for the calibration and measurement operations. When the user left-clicks the Float/Fixed button  on the Measurement Toolbar , the Measurement Toolbar will be fixed. In this case, the Camera Control Panel will not pop up automatically even if user moves mouse to the left side of the video windows. Only when the user left-clicks the  button on the Measurement Toolbar to exit from the measuring procedure will he be able to do other operations on Camera Control Panel , Autofocus Control Panel , or Synthesis Camera Control Toolbar . During the measuring process, when a specific measuring object is selected, an Object Location & Attributes Control Bar  will appear for changing location and properties of the selected objects.

6.2 The Camera Control Panel on the Left Side of the Video Window


Camera Control Panel	Function	Function Description
	Snap	Capture or Snap image from the current video window
	Record	Record video from the current video window
	Auto Exposure	Checking Automatic Exposure box will automatically adjust exposure time according to the Exposure Compensation value
	Exposure Compensation	Available when Auto Exposure is checked. Slide to left or right to adjust Exposure Compensation according to current video brightness to achieve proper video brightness
	Exposure Time	Available when Auto Exposure is unchecked. Slide to left or right to decrease or increase the exposure time to adjust the video brightness
	Gain	Adjust the gain value to decrease or increase the video brightness. The noise will be reduced or increased accordingly
	Red	Slide to left or right to decrease or increase the proportion of Red in the video window
	Green	Green is a base for reference and cannot be adjusted
	Blue	Slide to left or right to decrease or increase the proportion of Blue for the video
	White Balance	Auto White Balance adjustment according to the video window
	Sharpness	Adjust Sharpness level of the video window
	Denoise	Adjust Denoise level of the video window
	Saturation	Adjust Saturation level of the video window
	Gamma	Adjust Gamma level of the video. Slide to the right to increase the gamma value and to the left to decrease the gamma value.
Contrast	Adjust Contrast level of the video. Slide to the right side to increase and to the left to decrease video contrast	
DC	For DC illumination, there will be no fluctuation under the light source so no need for compensating light flickering	
AC(50HZ)	Check AC(50HZ) to eliminate flickering “strap” caused by 50Hz illumination	
AC(60HZ)	Check AC(60HZ) to eliminate flickering “strap” caused by 60Hz illumination	
Default	Set all the settings in the Camera Control Panel to the default values.	




The **Camera Control Panel** controls the camera to achieve the best image quality according to the specific applications; It will pop up automatically when the mouse is moved to the left side of the video window (in measurement status, the **Camera Control Panel** will not pop up. Only when the measurement process is terminated will the **Camera Control Panel** pop up by moving mouse to the left side of the video window). Left-clicking  button to achieve **Display/ Auto Hide** switch of the **Camera Control Panel**;

























6.3 The Measurement Toolbar On The Upper Side Of The Video Window



6.3.1 Introduction to Measurement Toolbar


The **Measurement Toolbar** will pop up when moving the mouse to any place near the upper side of the video window. Here is the introduction of the various functions on the **Measurement Toolbar**:



Icon	Function	Icon	Function
	Float/ Fix switch of the Measurement Toolbar	<input checked="" type="checkbox"/> Visible	Define measuring object in Show up/ Hide mode
	Select the desired Measurement Unit		
	Choose the same Magnification as the digital microscope current Zoom Ratio to ensure accuracy of measurement result when measurement unit is not in Pixel unit		

	Object Select		Point
	Angle		Four-point method to measure the angle
	Arbitrary Line		Three-Point method to measure the spacing
	Three-Point method to measure vertical line		Parallel Line
	Horizontal Line		Vertical Line
	Rectangle		Center + Radius Circle
	Three-points Circle		Ellipse
	Annulus		Two Circles
	Three-points Two Circles		Arc
	Polygon		Curve
	Arrow		Scale Bar
	Make Calibration to determine the corresponding relation between magnification and resolution, this will establish the corresponding relationship between the measurement unit and the sensor pixel size. The monitor's size can be input to achieve the accurate value of the digital magnification. The Calibration needs to be done with the aid of a ruler with an accuracy of more than 1mm. The detailed Calibration process is described in Sec. 6.3.2		
	Export the measurement information to CSV file(*.csv)		
	Delete all the Measurement Objects		
	Setting		Exit from Current Measurement Mode
	When the measurement ends, left-click on a single measuring object and the Object Location & Properties Control Bar will show up. The icons on the control bar mean Move Left , Move Right , Move Up , Move Down , Color Adjustment , and Delete .		

Note:1) When the user left-clicks **Display/Hide** button  on the **Measurement Toolbar**, the **Measurement Toolbar** will be fixed. In this case, the **Camera Control Panel** will not pop up automatically even if moving the mouse cursor to the left side of the video window. Only when users left-click the  button on the **Measurement Toolbar** to exit from the measurement mode will they be able to perform other operations in the **Camera Control Panel**, the **Autofocus Control Panel**, or the **Synthesis Camera Control Toolbar**.


2) When a specific measuring object is selected during the measuring process, the **Object Location & Attributes Control Bar**  will appear for changing the object location and properties of the selected objects.

3) To ensure accuracy of the measurement, after the calibration is turned on, the camera will automatically reset, and then sets the normalization magnification to 20X, and adjusts the focus to the required standard object distance. If the “**Calibration Object**” on the stage is not clear on the monitor, you need to manually adjust the height of the bracket to the clearest position, which is the standard object distance. After the **Calibration** is completed, use the **Measurement Toolbar** to measure the 1mm physical distance on the ruler, which should display 1mm on the monitor.

4) Even if the **Calibration** has been completed, once the user needs to measure, but is not sure whether the camera is at the standard object distance position, it is always better to reset it first, adjust the stand height in the reset state to make the observation object clear, and ensure that the camera is at the standard object distance position before measurement.

6.3.2 Calibration Method

User needs to prepare an **Calibration Object** such as ruler before **Calibration**;

Move the mouse to the upper side of the video window, the **Measurement Toolbar** will appear. Clicking  **Calibration** on the **Measurement Toolbar** to start the calibration.

1)The **XFCAMView** will pop up a message box: “1. Camera resetting for calibration...”

2)After the reset is finished, a message box: “2. Please put the calibration object on the stage(if not), adjust the height of the stand until the calibration object is in focus, then click OK button; ” will pop up.

3)After clicking the OK Button, XFCAMView will pop up a Calibration dialog shown below:




Figure 6-2 A Dialog for Calibration

Magnification:	the Magnification edit box, can be set from 1 to 20 as user want, Click Apply button to confirm;
Actual Length:	the Actual Length of the Calibration object on the stage, the unit can be selected with the right drop-down list box. Read the hint on the Calbration dialog to get the correct Calibration result;
Monitor Size:	the Monitor Size in Inch for the magnification calculation of the object displayed on the monitor;
Pixel:	the length in Pixel of the Calibration Line on the monitor;
Resolution:	the resolution in Pixel/Meter unit which is arrived by Pixel/Actual Length ;
OK:	Click OK button to end the Calibration ;
Users can refer to the message: “3. Please adjust the length and position of the calibration line until it is aligned with the calibration object. Then write down the actual dimension of the calibration object in the actual length edit box, and click OK.” to get the correct calibration result.	

The default monitor size is 27.0 inches. Users can enter the practical **Monitor Size**.

6.4 Synthesis Camera Control Toolbar At The Bottom Of The Video Window



Icon	Function	Icon	Function
	Zoom In the Video Window		Zoom Out the Video Window
	Horizontal Flip		Vertical Flip
	Color/gray		Video Freeze
	Display Cross Line		Image Overlay
	Auto Focus Control Panel		LED Brightness Control
	Compare Image with the Current Video		Browse Images Videos
	Settings		Check the Version of XFCAMView

The Setting function is relatively more complicated than the other functions. Here is more info about it:

6.4.1 Setting>Network

6.4.1.1 Setting>Network>General




Figure 6-3 Comprehensive Network General Settings Page

Name	The current camera name recognized as the network name
-------------	--

6.4.1.2 Setting>Network>LAN




Figure 6-4 Comprehensive Network LAN Settings Page

DHCP	Dynamic host control protocol allows DHCP server to automatically assign IP information to the camera. Only in Sec 6.4 LAN networking this item should be checked, so that cameras can automatically get IP information from routers/switches to facilitate networking operation;
Unicast/Multicast	By default, unicast function is used. Only in Sec 6.4 networking environment, when the router/switch has multicast function, camera can switch to multicast mode, which can save the network bandwidth consumed by the camera and facilitate the connection of more cameras in the same network;
IP Address	<p>Every machine on a network has a unique identifier. Just as you would address a letter to send in the mail, computers use the unique identifier to send data to specific computers on a network. Most networks today, including all computers on the Internet, use the TCP/IP protocol as the standard for how to communicate on the network. In the TCP/IP protocol, the unique identifier for a computer is called IP address.</p> <p>There are two standards for IP address: IP Version 4 (IPv4) and IP Version 6 (IPv6). All computers with IP addresses have an IPv4 address, and many are starting to use the new IPv6 address system as well.</p> <p>Users must manually configure their IP addresses on the camera side and computer side. The IP addresses set on the camera</p>

	side and computer side should be in the same network segment. The specific settings are shown Figure 6-. It's usually a private address. Private address is a non-registered address used exclusively within an organization. The internal private addresses retained are listed below: Class A 10.0.0-10.255.255; Class B 172.16.0-172.31.255.255; Class C 192.168.0-192.168.255.255. The suggested IP address is Class C.
Subnet Mask	Subnet Mask is used to distinguish network domain from host domain in 32-bit IP address;
Default Gateway	A default gateway allows computers on a network to communicate with computers on another network. Without it, the network is isolated from the outside. Basically, computers send data that is bound for other networks (one that does not belong to its local IP range) through the default gateway; Network administrators configure the computer's routing capability with an IP range's starting address as the default gateway and point all clients to that IP address.

Uncheck the **DHCP** and select the **Unicast** item, user still need to set the **IP address**, **Subnet** mask and **Default Gateway** as shown below:




Figure 6-5 Manual DHCP and Unicast

Uncheck the **DHCP** and select the **Multicast** item, user still need to set the **IP address**, **Subnet Mask** and **Default Gateway** as shown below:

6.4.1.3 Setting>Network>WLAN

Wi-Fi Mode	AP/STA mode to select;
Channel/SSID	Channel for the AP mode and SSID for the STA mode. Here, the SSID is the router's SSID ;
Password	Camera Password for the AP mode. Router Password for the STA mode




Figure 6-6 Network Setup

6.4.2 Setting>Measurement

This page is used for the define of the **Measurement Object** properties.





Figure 6-7 The Measurement Setup

Global	Used for setting digits behind the decimal point for measurement results;	
Calibration	Line Width	Used for defining width of the lines for calibration;
	Color	Used for defining color of the lines for calibration;
	EndPoint	Type: Used for defining shape of the endpoints of lines for calibration: Null means no EndPoint , rectangle means rectangle type of endpoints. It makes alignment more easily;
Point, Angle, Line, Horizontal Line, Vertical Line, Rectangle, Circle, Ellipse, Annulus, Two Circles, Polygon, Curve		
	Left-click the  along with the Measurement command mentioned above will unfold the corresponding attribute settings to set the individual property of the Measurement Objects .	

6.4.3 Setting>Magnification

This page's items are formed by the **Measurement Toolbar**'s **Calibration** command.




Figure 6-8 Comprehensive Magnification Settings Page

Name	Names such as 4X,10X, 20X, are based on magnification of the Digital microscopes.
Resolution	Pixels per meter. Image device like microscopes have high Resolution value;
Clear All	Click the Clear All button will clear the calibrated magnifications;

Delete	Click Delete to delete the selected magnification;
Up	Select a row in the magnification ratio and click Up to move up the currently selected magnification ratio;
Down	Select a row in the magnification ratio and click Down to move down the currently selected magnification ratio;

6.4.4 Settings>Image Format

Image Format	<p>JPEG: The extension of JPEG file can get very high compression rate and display very rich and vivid images by removing redundant images and color data. In other words, it can get better image quality with the least disk space. If measurement objects are available, the measurement objects will be burned into the image and the measurement cannot be edited.</p> <p>TIFF: TIFF is a flexible bitmap format mainly used to store images including photos and artistic images.</p>
Measurement Object Saving Method	<p>Burn in Mode: The measurement objects are merged into the current image. User could not edit the measurement objects any more. This mode is not reversible.</p> <p>Layered Mode: The measurement objects are saved in different layer with current image data in the target file. User could edit the measurement objects in the target file with some software on the PC. This mode is reversible.</p>




Figure 6-9 Comprehensive Image Format Settings Page

6.4.5 Setting>Video

Video Playback	Fast Forward/Reverse internal in second unite for Video Playback
Video Encode	Select the Video Encode format. Can be H264 or H265 . Compared with H264 , H265 has a higher H265 compression ratio which is primarily used to further reduce the design flow rate, in order to lower the cost of storage and transmission




Figure 6-10 Comprehensive Setting of Video page

6.4.6 Setting>Storage

<p>File System Format of the Storage Device</p>	<p>List the file system format of the current storage device</p> <p>FAT32: The file system of SD Card is FAT32. The maximum video file size of single file in FAT32 file system is 4G Bytes;</p> <p>exFAT: The file system of SD Card is exFAT. The maximum video file size of single file in FAT32 file system is 16E Bytes;</p> <p>NTFS: The file system of SD Card is NTFS. The maximum video file size of single file is 2T Bytes.</p> <p>Unknown Status: SD Card not detected or the file system is not identified;</p>
<p>Note: For USB Flash Drive, USB 3.0 interface is preferred.</p>	



Figure 6-11 Comprehensive Setting of Storage Page

6.4.7 Setting>Files




Figure 6-12 Comprehensive Setting of Files Name

<p>Image/Video File Name</p>	<p>Provide Auto or Manual naming paradigm for Image or Video file;</p>
<p>Auto</p>	<p>With specified name as the Prefix and XFCAMView will add digital after the Prefix for the Image or Video file;</p>
<p>Manual</p>	<p>A file dialog will pop up to enter the Image or Video file name for the captured Image or Video.</p>

6.4.8 Setting>Time




Figure 6-13 Time Setting

Time	User can set Year , Month , Day , Hour , Minute and Second in this page.
-------------	--

6.4.9 Setting>Language



Figure 6-14 Comprehensive Setting of Language Selection Setting Page

English	Set language of the whole software into English;
Simplified Chinese	Set language of the whole software into Simplified Chinese;
Traditional Chinese	Set language of the whole software into Traditional Chinese;
Korean:	Set language of the whole software into Korean;
Thailand	Set language of the whole software into Thailand;
French	Set language of the whole software into French;
German	Set language of the whole software into German;
Japanese	Set language of the whole software into Japanese;
Italian	Set language of the whole software into Italian;
Russian	Set language of the whole software into Russian;

6.4.10 Setting>Miscellaneous




Figure 6-15 Comprehensive Miscellaneous Settings Page

Clarity Factor	Check this will show the Clarity Factor on the video window screen to tell if the camera is focused correctly or not;
Ruler	Select to display the ruler in the video window, otherwise not to display the ruler;
Measurement	Select to display the measurement toolbar in the video window, otherwise not to display the measurement toolbar;
Overlay	Select to support saving graphics overlay information in fusion mode, otherwise it will not support;
Grids	Select to support saving mesh information in fusion mode, otherwise not to support;
USB video output switch back to mouse operation	Select automatic restart or manual restart to switch from USB video output to mouse operation;
ROI Color	Choosing the ROI rectangle line color
Cursor	Choosing the Cursor size according to the screen resolution or personal preference
Auto Exposure	Define the maximum automatic exposure time;
Auto Exposure Region	Select the AE reference area;
Camera Parameters Import	Import the Camera Parameters from the SD Card or USB flash drive to use the previously exported Camera Parameters
Camera Parameters Export	Export the Camera Parameters to the SD Card or USB flash drive to use the previously exported Camera Parameters
Reset to factory defaults	Restore camera parameters to its factory status;

6.4.11 Auto Focus Control Panel on the Right Side of Video Window

	Zoom Slider	Move the Zoom Slider to change the Zoom Ratio , the value will be displayed below the slider. It can be edited to set the desired Zoom Ratio
	Zoom Button	There are 3 Zoom Buttons , users can set specific zoom ratio for the quick control
	Optical Magnification	Optical Magnification is the designed lens magnification
	Digital Magnification	Digital Magnification is the object length on the monitor divided by the actual object length
	Focus Slider	Move the Focus Slider to change the focus lens position; The focus lens position value will be displayed below the slider. It can be edited to set the desired focus lens position;
Manual Focus	With Manual Focus radio button is checked, users can move the Focus Slider to change the focus lens position to get a clear image. The position value of the focus lens below the slider can be set by the user	

Autofocus	With Autofocus radio button is checked, the system will automatically focus the object on the stage, the focus lens position value under the Focus Slider will be refreshed in real-time; When the ROI or Object state is changed, the camera will perform the Auto Focus operation automatically
One Push	Clicking One Push button will perform a Autofocus operation at a time
Reset	Click Reset button to reset the Zoom and Focus modules. After the process is finished, the Zoom is set to 20X normalized magnification, and the Focus is fixed at the standard object distance(195mm in this model), if the object(such as a ruler for Calibration) is not clear, adjust the stand bracket to move the object to the standard object distance. Note: (see Measurement Toolbar>Calibration items for details).

6.4.12 Focus Region On The Video Window





Figure 6-16 Focus region

The **Focus Region** is used for selecting the region of interest for **Auto Focus** operation. When user clicks the  button on the **Synthesis Camera Control Toolbar**, the **Focus Region** will pop up as well with the **Autofocus Control Panel**. Users can click any part of the video window to select the focus region for **Auto Focus** operation.

When users close the **Autofocus Control Panel**, the **Focus Region** will be closed automatically.